

Cannulae



Cannulae



- Cannulae group consists two hygienic hoses and three universal tips operates automatically by microswitch when one of the cannulae is taken from the holder.
- All component can be replaced without any technical service and tools. The suction hoses has two different inner diameters which are 10 mm and 17 mm. The high quality and hygienic tubes can be removed and disinfected or cleaned easily.
- Cannulae group also involves washable filter system that holds large particles and inlays or precious metal filling.
- The cannulae group can be installed to any cabinet, wall or edge of table. It has two outlet which are 90° and 180° for different types installation.



AirVenturi

- Corpus Air Venturi is a low-cost extremely reliable appliance, small in size and easy to install.
- It requires no maintenance.
- Its performance is sufficient for low operating head and flow.
- You can use the air venturi when high vacuum system is out of order. It can be used at examination room where weaker suction is sufficient.
- Corpus Air Venturi can be installed inside any dental unit.